

Please be sure this checklist is inside the envelope with the dried blood spot testing cards.

- Unique identifier (i.e. Name, Date of Birth)
- DBS cards – each in their own coin envelope (unsealed)
- Coin envelopes places in Bitran bag
- Dessicant in Bitran bag (1 large pack or 3 smaller packs)
- Humidity card in Bitran bag

Storage: Check-off which one applies

- One week or less at room temperature
- More than one week at 4° C (or below)

Required contact info of person submitting tests

Name: \_\_\_\_\_

Program/Location: \_\_\_\_\_

Phone #: \_\_\_\_\_



10-120-7003 (IND - Rev. - 02/22)

Please be sure this checklist is inside the envelope with the dried blood spot testing cards.

- Unique identifier (i.e. Name, Date of Birth)
- DBS cards – each in their own coin envelope (unsealed)
- Coin envelopes places in Bitran bag
- Dessicant in Bitran bag (1 large pack or 3 smaller packs)
- Humidity card in Bitran bag

Storage: Check-off which one applies

- One week or less at room temperature
- More than one week at 4° C (or below)

Required contact info of person submitting tests

Name: \_\_\_\_\_

Program/Location: \_\_\_\_\_

Phone #: \_\_\_\_\_



10-120-7003 (IND - Rev. - 02/22)